

# **MEDICHOICE®**

## **WASH**

### **NEUTRAL pH INSTRUMENT CLEANER**

**MEDICHOICE®** Wash Neutral pH Instrument Cleaner is designed for use in ultrasonic cleaning systems, cart washers and for manual cleaning. It is formulated to remove soils and contaminants from instruments without causing degradation to their appearance. This is part of **MEDICHOICE®** total instrument care system process that prepares instruments for complete disinfection and sterilization.

#### **ADVANTAGES OF MEDICHOICE® WASH**

- Neutral pH
- Contains a wide range of anionic and nonionic surfactants
- Most effective and safe way that maximizes instrument cleanliness and service life
- Safe to use on brass, copper, aluminum, anodized aluminum and stainless steel
- Contains no phosphates
- Non-corrosive
- Rinse aid additive
- Reduces cost of instrument repair and replacement

**NOTE:** Packaged in four 1-gallon containers per case, and 15-gallon containers.